



NilRead

User Guide

October 2020

990-12027, Rev. 2.0

Copyright

© 2012-2020 Hyland Software Canada ULC. All rights reserved.

© 2007 The Closure Library Authors. All rights reserved.

Part Number: NIL-READ-GLB

Document Number: 990-12027 Rev. 2.0

Trademarks

Microsoft, Internet Explorer, Surface, ClickOnce and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Adobe, the Adobe logo, Acrobat, and Reader are either registered trademarks or trademarks of Adobe Systems Incorporated in the United States and/or other countries.

Firefox is a registered trademark of the Mozilla Foundation.

Chrome is a trademark of Google Inc.

iPhone, iPod, iPod Touch, iPad and Safari are trademarks of Apple Inc., registered in the United States and other countries.

Android is a trademark of Google Inc.

The Skype name, associated trademarks and logos and the “S” logo are trademarks of Skype or related entities.

FHIR is a registered trademark of HL7.

All other marks are trademarks or registered trademarks of their respective owners.

This page intentionally left blank.

Copyright	2
Trademarks	3
Intended Use	23
About NilRead	23
Version number	24
Warnings and precautions	24
General usage	24
Use on mobile devices	26
Measurements	26
Intended use within the USA	27
Supported modalities	27
Hardware requirements	27
Desktop computers	28
Mobile devices	28
Device specifics	28
Desktop computer with a diagnostic monitor	29
Mobile devices	29
Touch gestures	29
Reading environment verification	29
Studies	30
Open a study	30
Use presentations	33
Select presentations	34

Approve a presentation	35
Create a draft presentation	35
Delete a presentation	36
Use the patient timeline	36
View study or report details	37
View a study	37
View a report	37
Select priors to include in the timeline	38
Retrieve prior studies	39
Filter the timeline using one-click	39
Filter the timeline using Plus One (+1)	40
Example scenario	40
Timeline display options	41
Image stitching	41
Hide images on the screen	42
Add notes to a study	42
Apply presets	44
Create a series with all images	45
View stereometric images	45
View CAD marks	46
Use implant masks to reduce brightness	46
View DBT slice position	47
View study annotations	47
Edit videos	49

Modify a selection box	53
Add and remove selection boxes	54
View thumbnail images	54
Change the position indicator	55
Preview the edited video	56
Save the edited video	56
Return to the image viewer	57
Track review status	57
Organize studies in review folders	57
Manage reviewed studies	57
Use image visualization tools	58
Scroll	60
Pan	60
Zoom	61
Rotate	61
1:1 display mode	61
Configure your environment to use the 1:1 display mode	62
Enable the 1:1 display mode	62
Apply the 1:1 display mode	63
Quadrant Zoom	63
Invert	64
Window Level	64
Gamma	64
Enhance	64

Smart Zoom	65
Zoom In on a Portion of an Image	65
Compare Images	65
Relate	66
Link	67
Annotations and measurements	67
Propagate annotations and measurements	68
Annotations	68
Linear Measurements	70
Area Measurements	71
Angle Measurements	75
Color	76
Manage	76
Measurements for wide field ophthalmic photography images	76
Key Image	77
Titles	79
Full Quality	80
Reset	80
Hanging Protocols	80
Study Layout	80
View	81
Rendering	81
Reference	81
Thickness	81

Clipper	81
Curved MPR	82
Segment	87
NM Map	88
Fusion Map	88
Fusion Blend	88
Cine	88
RT Template	89
First, Previous, Next, Last	89
Use the image viewing area	90
Arrange images	90
Apply a study layout to the image viewing area	90
Apply a view to a screen	90
Sample layouts and views	91
1x2 layout	91
Multiple viewports	92
Multiple monitors	93
Side-by-side series comparison (study layout)	93
3D view	94
PET-C fusion view	94
Slab view	95
Multi-monitor auto-aligned mammography hanging protocol	95
Non-symmetric 2x1-1/3 study layout with ECG view	96
View the side panel	96

Select a series	97
Select hanging protocols	98
Apply a hanging protocol	98
Create a new hanging protocol	98
Edit a hanging protocol	99
Add a stage to a hanging protocol	99
Use full screen view	100
Maximize a viewport	100
Use reference lines	100
View image details	100
2D Images	101
Slab Images	104
3D Images	107
View image orientation	109
View DICOM attributes	110
View study information	110
Use worklists and folders	110
About worklists and folders	110
Worklists	111
Folders	111
Manage worklists	111
Add a worklist	112
Edit or delete a worklist	115
Manage folders	116

Add a folder	116
Edit or delete a folder	118
Manage the studies in a folder	119
Add studies to a folder	119
View and edit folder notes	119
Organize studies by date	119
Move studies to another folder	120
Remove studies from a folder	120
Open a study in a worklist or folder	120
Manage studies	121
Access studies	121
Use Patient Search	122
Open multiple studies	124
Use “Break Glass” to find studies	124
Lock studies	125
Upload studies	125
Upload DIOCM files from media	125
Upload zipped DICOM files	126
Upload a folder containing DICOM files	126
Upload individual DICOM files	126
Upload non-DICOM files	127
Create an upload link	127
Attach files to studies or patients	128
Edit or split studies	129

Download studies, series or images	129
Merge patients	130
Merge studies	131
Delete studies, series or images	131
Delete series containing segmentation results	131
Export images	132
Print images	132
View reports	133
Review ECG studies	134
About ECG paper	134
Use ECG views	135
Understanding QRS Axis Interpretation	138
Use the ECG view tools	139
Gain and Timescale	139
Filter	141
HR	142
QTQTc	144
Annotation Tools	146
Change ECG preferences	147
Use the vessel tracer	147
About the vessel tracer	147
Use the Vessel Tracer view	147
Define a vessel	149
Adjust the centerline and wall	151

View radiation therapy (RT) plans	153
Using the RT Plan panel	153
Structures	159
Isodoses	163
Beams	175
About RT plans	177
Using the RT Graph view	178
Using the DRR view	182
Apply an RT template to a study	185
Manage RT templates	186
Add a template	186
Edit or delete a template	191
Enable a template	191
Import and export templates	192
Anonymize patient data	192
About anonymization	192
Manage confidentiality profiles	193
Add a confidentiality profile	193
Edit or delete a confidentiality profile	194
Refresh the confidentiality profiles list	195
Manage confidentiality masks	195
Add a confidentiality mask	195
Add regions to be redacted	196
Edit or delete a confidentiality mask or region	196

Refresh the confidentiality masks or regions list	196
Left Ventricle Ejection Fraction (LVEF)	197
Selecting the end diastole	198
Measuring the volume using the Simpson method	198
Selecting the end systole	199
Improving precision using biplane Simpson volume	200
Viewing all volumes used for ejection fraction calculation	200
Replacing one volume with another	201
Performing another ejection fraction measurement in the same study	201
Collaborate	202
Use online meetings	202
Create a meeting for a folder	203
Start a meeting for a folder	203
Send the meeting link to participants	203
Join a meeting for a folder	204
Create a general meeting	204
Start a meeting	204
Join a meeting	204
Use meeting controls	205
Use Skype	206
Send secure study links	207
Send a secure link by email	207
Copy a link	208
Send series links	208

Share bookmarked images	210
Create a bookmark	210
View a bookmark	210
Send a bookmark	211
Delete a bookmark	212
Share secondary capture images	212
Create a series of secondary capture images	212
Share a series of secondary capture images	213
Delete a series of secondary capture images	214
Preferences	214
Manage your user profile	214
Change your password	215
Manage study note templates	216
Add a template	216
Edit or delete a template	223
Enable a template	223
Import and export templates	223
Manage assets	223
Add an asset	224
Edit or delete an asset	224
Change your user preferences	224
Change your workstation preferences	226
Change mouse, keyboard and touch preferences	228
Step 1: Add a template	228